# User's Manual

## EUSRA RF Electrode

<b>EN</b>	<b>Electrode for electrosurgical device</b>	1 page
<b>DA</b>	<b>Elektrode til elektrokirurgisk indretning</b>	8 page
<b>DU</b>	<b>Elektrode voor elektrochirurgisch apparaat</b>	15 page
<b>FR</b>	<b>Électrode pour dispositif électrochirurgical</b>	22 page
<b>GE</b>	<b>Elektrode für elektrochirurgische Geräte</b>	30 page
<b>GR</b>	<b>Ηλεκτρόδιο για ηλεκτροχειρουργική συσκευή</b>	38 page
<b>IT</b>	<b>Elettrodo per dispositivo elettrochirurgico</b>	47 page
<b>NO</b>	<b>Elektrode for elektrokirurgisk enhet</b>	55 page
<b>PO</b>	<b>Elétrodo para dispositivo electrocirúrgico</b>	62 page
<b>SP</b>	<b>Electrodo para dispositivo electroquirúrgico</b>	70 page
<b>SW</b>	<b>Elektrod för elektrokirurgisk anordning</b>	78 page
<b>TU</b>	<b>Elektrocerrahi cihazı için elektrot</b>	85 page

### User's Manual

#### 1. Introduction

The electrode for this electrosurgical device is used by connecting it to the RF lesion generator from STARmed Co., Ltd. which was designed to coagulate tissues.

The electrode, in combination with an endoscopic ultrasound (EUS) or scope, is used to access body cavities in order to coagulate tissues. Tissues can be coagulated with the supply of radiofrequency energy from the electric cautery, which is connected to the generator for monitoring the impedance and temperature of tissues.

The electrode's fixed insertion length can be adjusted by the screw on the handle from 0 cm to 8 cm (1 cm units).

The 19G Model electrode length can be adjusted from 1 cm to 5 cm by the controller which adjusts the length of endoscope.

#### 2. Guide for Use

Read all of the following information and also the provided user's generator guide before use.

The EUSRA RF Electrode is designed to be used in conjunction with the RF lesion generator from STARmed Co., Ltd. If the guidelines are not properly followed, damages can be caused from electricity or heat and the equipment may malfunction.

#### 3. Intended Purpose

Sterile single-use RF coagulation electrode is intended for coagulation of tissue during endoscopic, percutaneous, laparoscopic, and intraoperative surgical procedures.

#### 4. Contraindications

Sterile single-use RF coagulation electrode is

contraindicated for use in patients with Implantable Pacemakers and Automatic Implantable Cardioverter/Defibrillators (AICDs), as they may be adversely affected by radiofrequency (RF) current.

#### 5. Intended patient population

age	Newborn to geriatric
weight	>2.5Kg
gender	Male and Female
health	Use for patients with liver cancer, thyroid cancer, lung cancer, kidney cancer or other lesions.
Nationality	Multiple
Patient state	The patient is not the device user. Not relevant unless the patient is agitated.
Contraindication group	Patients with implantable pacemaker and implantable cardioverter/defibrillators.

#### WARRANTY

STARmed Co., Ltd. warrants sufficient care for the design and manufacture of this device. We do not warrant any content that is not described in this manual. The handling, storage, and operation of the EUSRA RF Electrode beyond what is described in the IFU may directly affect the results obtained from its use.

STARmed's device under this warranty is limited to repair and replacement, and STARmed Co., Ltd is not responsible for any economic expenditures caused by unexpected issues, significant loss, or damage. STARmed Co., Ltd. is not transferring responsibility and duty related to this device to other parties, and will not take any responsibility regarding reuse and usage of re-sterilized or expired products.

#### 6. Warning

- Before using this product, check whether there is any damage to the packaging.

- The expiration date is three (3) years from the sterilization date. Use of the products after the expiration date is prohibited.
- This device is disposable and reuse is prohibited.
- This device is designed to be used by a qualified and trained medical doctor.
- This device should be stored in a cool place without sunlight and humidity.
- Do not touch the electrode's active tip while the power supply is on.
- Do not touch the shaft of the electrode with fingers or other tools (insulator and components are made of metallic).
- Do not touch the electrode's hand piece/cable with metallic materials. Touching the cable with metallic materials may cause the patient or operator to get shocked, cause a fire, or have an accident with injury.
- Do not use on patients with cardiac pacemakers or other active implants
- Do not use near conductive materials.
- Like all other electrosurgical devices, do not use in the presence of flammable anesthetics, nitrous oxide or any other flammable substances. Doing this may cause a fire in the electrosurgical device.
- Flammable gases or other materials with a high risk of flammability should be prohibited from being placed nearby the RF generator. If a flammable detergent has been used for cleaning or disinfection, the liquid should be evaporated before the start of the procedure.
- Cloth and gauze soaked in oxygen can possibly be ignited from sparks which occur from normal use of the device. Please use with caution.
- Do not activate the electrode while it is touching any metallic substances or tools. Doing this may cause an unexpected injury to the patient or operator.

- Do not touch metals or tools with the electrode while the power supply is on. Doing this may cause an unexpected injury to the patient or operator and/or damage to the electrode or other devices.
- When using the electrode, do not hold the tip of the electrode by forceps, do not pass the electrode through a hole of other devices, and do not scratch the surface on the electrode's tip by sharp objects. These actions might destroy the insulation of the electrode and/or damage the electrode. It might cause unexpected injury to the patient or the user.
- Do not access the lesion without proper imaging guidance. It might cause an unexpected injury to the patient and/or damage to the electrode.

### 7. Potential Complications

- Burning by overheating of the surgical device
- Dangers from an inexperienced operator's use
- Side effects or cross infection from reuse
- Weakness of liver functions
- Delayed bleeding in the operated body parts
- Recurrence of cancer
- Perforation
- Respiratory depression or arrest
- Symptoms after RFA treatment including: abdominal pain, fever, nausea, headache, right shoulder joint pain, and chest discomfort

The potential complications are not limited to the list above. Other complications may occur.

### 8. Before Use

(1) Preparation of coolant

- The temperature of the coolant should be maintained at almost 0°C by keeping the IV bag in the refrigerator for four (4) or more hours before the procedure. It should be used just before performing the procedure.

**⚠ Caution:** Using coolant from an external and/or openable container may cause the electrode to be clogged by an inflow of foreign substances.

### (2) Grounding pads

**⚠ Caution:** Re-using grounding pads is prohibited.

**⚠ Caution:** Attach an even number of grounding pads to the patient. If not, burns may occur from the grounding pad(s).

- Attach the horizontal lines of the grounding pads on the thighs widthwise after removing each grounding pad's transparent plastic protection film.

- Be careful not to form air bubbles between the grounding pads and the skin while attaching. If required, remove hair at the site of treatment. Clean and dry the area before attaching the grounding pads.

- Attach the grounding pads at the same distance from the treated area.

- Make sure the lines connected to the grounding pads and the electric surgical device are not twisted.

- Recommended attachment site of the grounding pads: muscular body area with large heat dissipating blood vessels

**⚠ Caution:** The areas to be avoided for attachment of the grounding pads: areas with metal implants, a pacemaker, fat, previous burns, inflammation, bone protrusion(s), an ECG electrode and/or electric line(s), and/or areas where liquid can be collected.

### (3) Tubing set-up

- Connect the electrode, the inflow/outflow tubes, and the pump tubing.

**⚠ Caution:** Be careful when taking out the electrode from the packaging. The product's elasticity may cause injury to the patient and/or operator.

- Open the cover of pump head upwards, place the pump tube (rubber tube) in the middle of the roller, and place the tube precisely at the groove of the pump head. Then, secure the tube in the pump by closing the cover.



- Straighten the tube by stretching both sides.

- Hang the coolant on the IV pole and puncture the bag with the spike. The roller clamp should be temporarily closed until the beginning of the procedure.

- Open the roller clamp during the procedure.

### (4) Length Control Luer ( Only to 18G Model)

- Use the Length Control Luer to change the length depending on the length of the EUS or scope.

## 9. During Use

### 9-1. 18G model

(1) Connect the electrode to the RF lesion generator. Confirm the connection status

**⚠ Caution:** Before inserting the electrode into the

lesion, check if the exposed length of the electrode and other specifications are identical with the specifications on the labels.

(2) Check the location of the EUS or scope when targeting the lesion for accuracy.

 Caution: Refer to the user's manual for warnings and how to use the EUS or scope.

(3) The fixed part of the insertion length should be locked by the screw at the point of zero (0) before the electrode is inserted into the instrument channel of the EUS or scope.

 Caution: Keep the electrode tip in the sheath when combining it with the EUS or scope.

(4) Connect the connection luer with the instrument channel luer by turning the handle of the electrode after the electrode is put in the instrument channel of the EUS or scope.

 Caution: Do not change the tip length before the electrode is connected to the EUS or scope. Doing this may cause damage to the EUS or scope.

(5) Adjust the insertion length of the electrode by loosening the screw on the insertion length controller and moving the controller onto the target position. Then, tightly lock the screw on the controller.

 Caution: Inserting the electrode's tip without tightly locking the handle screw on the insertion length controller may cause the electrode tip to be inserted into an undesired tissue location. This may lead to perforation or injury by an undesired moving of the insertion length controller.

(6) Check the location of the lesion to be treated using the EUS or scope, and place the electrode in the lesion.

(7) Once the electrode is placed in the lesion, check that the coolant is flowing out from the electrode's output tube by operating the pump connected with the input tube. Turn on the RF lesion generator by

pressing the START button.

 Caution: Do not start the radiofrequency output before the inside of the electrode tip gets cold. This may cause carbonization.

(8) Once the procedure is completed, turn off the POWER switch of the RF lesion generator and pump.

(9) Pull back the tip of the electrode into the sheath and tighten the screw at the point of zero (0). The electrode can now be removed from the instrument channel of the EUS or scope.

(10) Be careful to remove the disconnected electrode from the instrument channel.

 Caution: Be careful that the patient and the operator are not injured by the electrode's tip when removing the electrode from the EUS or scope.

(11) Follow the user's manual when removing the EUS or scope from the patient.

### 9-2. 19G model

(1) Connect the electrode to the surgical instrument and then check whether it is properly connected.

 Caution: Before inserting the electrode into the lesion, make sure whether the exposure length of the electrode and other specifications match the label.

(2) Insert the EUS and Scope into the lesion to be treated, and then check the location of the lesion.

 Caution: For instructions and precautions about the use of EUS and Scope, see the user's manual of the corresponding products.

(3) Before inserting the electrode into the instrument channel of the EUS and Scope, check whether the insertion length locking part is placed at "0".

 Caution: When the electrode is inserted into the EUS and Scope while the tip of the electrode comes out of the sheath, it may damage the EUS and Scope.

(4) Put this product into the ultrasound endoscope in order to determine the appropriate length of the outer jacket depending on the endoscope.

(5) Insert the lure connected to the endoscope into the entrance of the endoscope working channel and then turn it.

(6) Lift the button of the controller which adjusts the length of endoscope and turn it 90 degrees to be moved.

(7) After that, adjust it until it reaches the appropriate length. Turn the button 90 degrees to be fixed.

 Caution: If the needle length of the electrode is increased or adjusted while the connection between the electrode and EUS and Scope is not ensured, the EUS and Scope may be damaged.

(8) Lift the button of the insertion length locking part and turn it 90 degrees to be moved.

(9) Move the button to the desired insertion length of the tip and turn it 90 degrees to be fixed.

 Caution: When the electrode TIP is inserted by advancing the electrode handle without fixing the button of the insertion length locking part, the insertion length locking device may move, and then the electrode TIP may be inserted into the tissue of the undesired portion, which may cause perforation or injury.

(10) Check the location of the lesion to be treated by using the EUS and Scope, and then advance the electrode handle to place the electrode TIP in the lesion.

(11) When the electrode is placed in the lesion, operate the pump connected to the inflow tubing to check whether the cold liquid comes out from the outflow tubing of the electrode. After that, press the power switch of the electrosurgical instrument to be

operated.

 Caution: If the high frequency is produced while the inside of the electrode TIP is not sufficiently cooled down, carbonization may take place in the electrode TIP, which may cause problems in product performance.

(12) When the operation is completed, turn off the power switches of the electrosurgical instrument and the pump.

(13) Pull the electrode handle backwards, and put the electrode TIP in the sheath. Set the insertion length locking part to "0", and fix the button of the insertion length locking part to the appropriate position. After that, turn the electrode handle to separate it from the instrument channel of the EUS and Scope.

(14) Carefully remove the separated electrode from the instrument channel.

 Caution: When the electrode is removed from the EUS and Scope, be very careful to avoid any damage in the patient and operator due to the electrode tip.

(15) According to the instructions of the product, remove the EUS and Scope from the patient.

### 10. After Use

(1) Grounding pad

- When the procedure is finished, turn off the POWER switch of the RF lesion generator and disconnect the grounding pad connectors from the generator.

- Slowly remove the grounding pads in order to prevent injury on the surface of the skin.

- Do not remove the grounding pads by pulling the electric cables of the ground pads.

(2) Turn off the POWER switch of the RF lesion generator and remove the electrode and other parts after the procedure is finished.

(3) Dispose of the electrode according to the waste materials disposal procedure.

### 11. Storage

(1) Store at room temperature.

 Caution: This device is intended to be sold and used by the order of a medical doctor according to the related medical device laws.

### 12. Action for Product Damage

 Caution: If any damage on the product is visually noticed, do not use to prevent injury to the patient and/or operator.

**13. Expiration Date:** Three (3) years from sterilization

### 14. Symbol

Authorized Representative in the European Community

Consult instructions for use

Manufacturer

Catalogue number, Model name

Sterilized using ethylene oxide

Caution

Use-by date

Batch code

Do not re-use

Date of manufacture

### Reuse Precaution Statement

The contents of this product have been sterilized with E.O. gas. If the sterilized package is damaged, do not use it and contact STARmed Co., Ltd.

This product can be used only one time. Reuse, re-treatment, or re-sterilization is not allowed.

If the electrode is reused, re-treated, or re-sterilized, it may cause defects in the structural functions of the electrode and damage of the RF lesion generator which may lead to injury, disease, or death to the patient. Reuse, re-treatment, or re-sterilization can cause device contamination and infections among patients or from other diseases. Contamination of this device may be a cause of injury, disease, or death.

### 15. Descriptions

Common name: Electrode for electrosurgical device

Trade name: EUSRA RF Electrode

Catalogue number, Model name: Refer to the product

Lot number: Refer to the product label

Manufactured date: Corresponding year and date

Number of items: 1 set

Sterilization method: E.O. gas

Expiration date: Three (3) years from date of sterilization

Method of use: Refer to the manual

Intended use: Refer to the manual

Cautions in use: Refer to the manual

Method of storage: Store at room temperature

Manufacturer: STARmed Co., Ltd.

**※ This product is a disposable medical device and reuse is not allowed.**

STARmed Co., Ltd.

Address: (Jungsandong, Daebang-Triplaon Business Tower), B-dong, 4F & 12F, 158, Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea

T E L: +82 506 816-3546

F A X: +82 506 816-4546

<http://www.STARmed4u.com>

E-mail: [inquiry@STARmed4u.com](mailto:inquiry@STARmed4u.com)

Electrode for electrosurgical device  
[EUSRA RF Electrode]

CE  
1639

**EC REP** EC Authorised Representative

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Üretici firma: STARmed Co, Ltd
STARmed Co, Ltd