# Niti-S Hot Nagi™ Stent & Electrocautery Stent delivery system

## Table of Contents

|                  |    |
|------------------|----|
| English.....     | 3  |
| Français.....    | 8  |
| Deutsch .....    | 14 |
| Italiano .....   | 20 |
| Español .....    | 26 |
| Nederlands ..... | 32 |
| Português .....  | 38 |
| Türkçe.....      | 44 |
| Čeština .....    | 49 |
| Slovensky.....   | 55 |
| Dansk.....       | 61 |
| Ελληνικά.....    | 67 |
| Lietuvių.....    | 73 |
| Polski .....     | 79 |
| Română.....      | 85 |
| Svenska.....     | 91 |
| Norsk .....      | 96 |

**MRI Information**

**MR Conditional**

The Niti-S Hot Nagi™ Stent is MR conditional. Non-clinical testing demonstrated that the Niti-S Hot Nagi™ Stent is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla or 3-Tesla, only
- Maximum spatial gradient magnetic field of 5,000-Gauss/cm (extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence)
- Normal Operating Mode of operation for the MR system

**MRI-Related Heating**

In non-clinical testing, the Niti-S Hot Nagi™ Stent produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

|                                                      | <u>1.5-Tesla</u> | <u>3-Tesla</u> |
|------------------------------------------------------|------------------|----------------|
| MR system reported, whole body averaged SAR          | 2.9-W/kg         | 2.9-W/kg       |
| Calorimetry measured values, whole body averaged SAR | 2.1-W/kg         | 2.7-W/kg       |
| Highest temperature change                           | 2.9°C            | 2.4°C          |
| Scaled to whole body averaged SAR of 2.0-W/kg        | 2.0°C            | 1.7°C          |

**Artifact Information**

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Niti-S Hot Nagi™ Stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of this implant. The lumen could not be visualized using these pulse sequences.

|                   |          |               |          |               |
|-------------------|----------|---------------|----------|---------------|
| Pulse Sequence    | T1-SE    | T1-SE         | GRE      | GRE           |
| Signal Void Size  | 702-mm2  | 317-mm2       | 883-mm2  | 482-mm2       |
| Plane Orientation | Parallel | Perpendicular | Parallel | Perpendicular |

# English

## Instructions for Use

### 1. Description

The Niti-S Hot Nagi™ Stent consists of an implantable metallic stent and an Electrocautery Stent Delivery System [HDS-10].

The stent is made of Nitinol wire. It is a flexible, fine mesh tubular prosthesis which has radiopaque markers on each end and at the center.

The Electrocautery Stent Delivery System is designed to penetrate the adjacent tissue of the human organs with electrocautery connected to the Electro Surgical Unit.

| Model Name             |
|------------------------|
| Niti-S Hot Nagi™ Stent |

Figure 1. Stent Model

The Stent is loaded in the Electrocautery Stent Delivery System [HDS-10] and upon deployment the body of the stent imparts an outward radial force while both flares impart a horizontal inward force allowing transmural drainage.

Niti-S Hot Nagi™ Stent is Fully Silicone Covered and can be removed (see Section 5, Warnings and Precautions).

The Electrocautery Stent Delivery System [HDS-10] is an electrocautery-enhanced delivery system that is compatible with therapeutic echoendoscopes having a working channel of **3.7 mm diameter or larger**.


Figure 2. Electrocautery Stent Delivery System [HDS-10]

- The endoscopic Electrocautery Stent Delivery System [HDS-10] has a usable length of 180cm

**Endoscopic** Type is recommended

- When approached endoscopically

### 2. Principle of Operation

The outer sheath is pulled back by immobilizing the active connector in one hand, grasping the active handle with the other hand, and gently sliding the active handle along the pusher towards the active connector. Retraction of the outer sheath releases the stent.

### 3. Indication for Use

The Niti-S Hot Nagi™ Stent and Electrocautery Stent Delivery System [HDS-10] is intended for the drainage of a pancreatic pseudocyst through a transgastric or transduodenal approach.

### WARRANTY

Taewoong Medical Co., Ltd. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Taewoong's control directly affect the instrument and the results obtained from its use. Taewoong's obligation under this warranty is limited to the replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

### 4. Contraindication

This device is contraindicated for use in any and all cardiovascular applications.

Additional contraindications include:

- Hemodynamic instability
- All others than indication for use
- Recapturing a stent during its deployment is contraindicated.
- Cystic neoplasms
- Pseudoaneurysms
- Duplication cysts
- Non-inflammatory fluid collections
- Patients with abnormal coagulation or who require ongoing complete anticoagulation at the time of implantation and post stent placement due to increased possibility of bleeding
- Patients with altered anatomy that precludes the physician's ability to deliver the stent. (e.g. too small lesion for the stent placement)
- Patients that have allergies or are sensitive to any of the device material
- Patients with contraindications to use electrical devices

## 5. Warnings and Precautions

Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- The safety and effectiveness of this device for use in the vascular system has not been established.
- The device should be used with caution and only after careful consideration in patients with radiation colitis or proctitis.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the Electrocautery Stent Delivery System [HDS-10] to organic solvents (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- Long-term patency of this stent has not been established and regular observation is recommended.
- The stent cannot be resheathed once deployment has been initiated.
- Do not attempt to recapture/reload a stent once its deployment is advanced.
- Care should be taken when removing the Electrocautery Stent Delivery System [HDS-10] and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after stent deployment, as this may result in perforation, bleeding, stent dislodgement or stent migration.
- The packaging and the device should be inspected prior to use.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The Niti-S Hot Nagi™ Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The Niti-S Hot Nagi™ Stent is intended for single use only. Do not resterilize and/or reuse the device.
- In case any defect is found, such as a foreign substance, stop using immediately and contact Taewoong Medical.
- Before use, examine the outer surface of devices which is intended to be inserted into a patient or used during procedure. Do not use a device that has unintended rough surfaces, sharp edges or protrusions which may cause harm. Cut, burned or damaged device insulation may cause unsafe current in either patient or operator.
- The maximum rated input voltage for this device is 1kVp-p.
- Select an appropriate EMC compatible electrosurgical unit (ESU).
- Any electrosurgical accessory constitutes a potential electrical hazard to the patient and operator.
- Before using this device, follow recommendations provided by the manufacturer of the electrosurgical unit for the proper placement and utilization of the patient return electrode. Ensure that a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.
- Examine all components before use. Do not use a device that has been cut, burned or damaged. Damaged device insulation may cause unsafe currents in either patients or operator.
- Power off the electrosurgical unit when not in use.
- When applying current, ensure the Electrocautery Stent Delivery System [HDS-10]'s active tip is completely outside the endoscope. Contact between the active element and the echoendoscope may cause grounding, which can result in injury to patient and operator, or damage to the endoscope.
- Do not touch the Electrocautery Stent Delivery System [HDS-10] pusher when applying current. It may cause electrical hazard to the patient and operator.
- Do not use with wet hands.
- Electrocautery Stent Delivery System [HDS-10]'s contents supplied STERILE using ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call manufacturer.
- Make sure inner tube of the delivery system is resheathed into the outer tube when remove the device after use.
- Unintended overtemperature around 48.2°C was measured during the Electrocautery Stent Delivery System. Do not use cautery for over 40 seconds to ensure patient safety.
- Do not use other than as stated in the instructions for this device and designated accessories.
- Do not use this device in any echoendoscope with a working channel smaller than 3.7mm.
- No modification of this equipment is allowed.
- Do not use in patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional (e.g. cardiologist). A possible hazard exists since interference with the electronic implant may occur, or the implant may be damaged.
- Interference of high frequency medical electrical equipment may adversely influence operation of other electronic equipment.
- Before use, compatibility with electrosurgical generators, accessories and other endoscopic equipment should be checked for safe use. Using incompatible equipment or equipment not specified by this Instruction for Use can result in patient injury or equipment damage.
- Do not use in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N<sub>2</sub>O) and oxygen) or in close proximity to volatile solvents (such as ether or alcohol), as explosion may occur.

- Do not place instrument near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause fire.
- Do not activate the instrument when it is not in contact with target tissue, as this may cause injuries due to capacitive coupling.
- The surface of the active electrode may remain hot enough to cause burns after the RF current is turned off.
- Use pure cut generator settings with Electrocautery Stent Delivery System [HDS-10]. Do not use blended, or coagulation generator modes. Blended or coagulation modes may result in failure to access the cyst, prolonged time to access the cyst, tissue tenting or resistance.
- Ensure correct generator installation. The generator must be installed and put into service according to the EMC information provided in the generator manufacturer's guidance and declaration for electromagnetic compatibility.
- Connect the Electrocautery Stent Delivery System [HDS-10] to the electrosurgical unit only when the unit is off. Failure to do so may result in an injury or electrical shock to the patient or operating room personnel.
- Due to concerns about the carcinogenic and infectious potential of electrosurgical by products (such as tissue smoke plume and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used during the procedure.
- Prior to increasing the ESU intensity, check the adherence of the return electrode and its connections. Apparent low output or failure of the device to function correctly at the normal operating settings may indicate faulty application of the return electrode or poor contact in its connections.
- Device must be used in conjunction with a Type BF or CF generator, see compatible electrosurgical unit or generator information.

## 6. Instructions in the Event of Damage

**WARNING:** Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

## 7. Potential Complications

Potential complications associated with the use of Niti-S Hot Nagi Stent may include, but are not limited to:

- Inadequate expansion
- Pancreatitis
- Abscess formation
- Haemorrhage
- Vomiting
- Peritonitis
- Hematoma
- Rupture of intracystic artery
- Stent removal failure
- Anesthesia complications
- Improper placement ; incomplete deployment; stent migration into the target structure or, GI Tract ; separation of coating material from stent; stent fracture ; coating material wear; coating material failure; puncture of coating material
- Tissue ingrowth or overgrowth leading to difficulty or a failure to remove stent
- Stent dislodgement
- Adverse reaction to implant materials and/or delivery system (e.g., abdominal or back pain, nausea, infection, fever, chronic inflammation or foreign body reaction)
- Minor or excessive bleeding requiring intervention
- Leakage of pseudocyst or bowel or biliary tract contents causing inflammation or peritonitis
- Stent occlusion
- Local infection at the implant site.
- Tissue damage during stent implantation and/or removal
- Ulceration or erosion of mucosal or organ wall linings
- Pneumoperitoneum
- Sepsis (bacterial, endotoxin or fungal)
- Perforation
- Surgical intervention(endoscopy, transfusion or surgery)
- Persistent connection to the target structure after removal (fistula)
- Unintended electrical shock, muscle stimulation or burns
- Cardiac arrhythmia or arrest
- Death

## 8. Application

- Environment
  - √ General: Operating room in the hospital
  - √ Physical condition: The user should maintain the conditions as below during the procedure.

|                       | Operation         | Transportation    |
|-----------------------|-------------------|-------------------|
| Range of temperature  | 10 °C to 40 °C    | -18 °C to 40 °C   |
| Range of humidity     | 10 % to 75 %      | 10 % to 90 %      |
| Range of air pressure | 70 kPa to 106 kPa | 70 kPa to 106 kPa |

## 9. Technical Specifications

|                                                     |                                                                                                                                    |
|-----------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| <b>Use</b>                                          | Sterile, Single-Patient Use                                                                                                        |
| <b>Energy</b>                                       | Monopolar                                                                                                                          |
| <b>Maximum rated Input</b>                          | 1.0kV peek-to-peek (500Vp)                                                                                                         |
| <b>Recommended Generator Settings</b>               | Pure cut mode, 80-120 Watts (400-500Vp)                                                                                            |
| <b>Compatible Electrosurgical Unit or Generator</b> | ERBE VIO Model                                                                                                                     |
| <b>Electrosurgical Generator Requirements</b>       | *Rated supply voltage: 100V - 120V / 220V - 240V<br>*Rated supply frequency: 50 / 60Hz<br>*Line current: 8A / 4A<br>*Pure cut mode |
| <b>Connectors</b>                                   | Monopolar endoscopic cable, 3mm female plug, Select cables specified by generator manufacturer.                                    |
| <b>Dispersive Pad</b>                               | Select pad or return electrode specified by generator manufacturer.                                                                |
| <b>Guidewire Compatibility</b>                      | 0.035 inches (0.89mm) insulated                                                                                                    |
| <b>Echoendoscope Compatibility</b>                  | Working channel of 3.7mm diameter or larger                                                                                        |

## 10. Preparation

Ensure that a 0.035 inch guidewire has been placed through the echoendoscope working channel and into the target structure. The echoendoscope elevator should be in the lowered (open) position.

## 11. Procedure

### - (180cm Electrocautery Stent Delivery System [HDS-10] usable length)

Endoscopy Ultrasound (EUS) should be performed prior to placement of the Niti-S Hot Nagi™ Stent to characterize the extent of the lesion and its morphology.

### ① Stent Size Determination

- The stent length and diameter should be determined by the doctor after endoscopic and/or fluoroscopic inspection of the lesion.
- The length and diameter of the stent should be chosen in a way that both walls of the created transmural drainage remain tightly together, thus preventing any migration.

### ② Stent Deployment Preparation

- After removing the stylet, slowly insert the Electrocautery Stent Delivery System [HDS-10] into the working channel of the echoendoscope over the guidewire. Connect the active cable(monopolar) to electrosurgical generator. (Do not inject saline solution into the active connector.)
- With the generator powered off, connect plug of the Electrocautery Stent Delivery System [HDS-10] to the active cable(monopolar). Power on the generator and ensure the correct settings (see recommended settings in Technical Specification Table). ONEY pure cut settings should be used (**80-120 Watts, 400-500Vp**). Using the ERBE electrosurgical generators, the “**Auto cut**” setting corresponds to the pure cut setting. Do not use ERBE “Endo cut” setting.
- Position the Electrocautery Stent Delivery System [HDS-10] to the target structure.
- Energize and advance. Verify the pure cut mode setting, energize the device and advance carefully into the target structure.
- After the device has entered into the structure, power off the generator and unplug the Electrocautery Stent Delivery System [HDS-10]



Figure 3

### ③ Stent Deployment Procedure

**PRECAUTION:** Do not twist the Electrocautery Stent Delivery System [HDS-10] or employ a boring motion during the deployment as this may affect positioning and ultimate function of stent.

- Under the fluoroscope and/or endoscopic guidance, position the Electrocautery Stent Delivery System [HDS-10]
- Once the Electrocautery Stent Delivery System [HDS-10] is in the correct position for deployment, unlock the proximal valve of the active handle by turning the valve more than twice in an anti-clockwise direction.
- To begin stent deployment, immobilize the the active connector in one hand and grasp the active handle with the other hand. Gently slide the active handle back along the pusher towards the active connector

- d) Under EUS and/or fluoroscopic guidance, deploy the distal flare of the stent, and then deploy the proximal flare of the stent under endoscopic guidance while making sure that the stent connects both walls together.

**Figure 4**

**CAUTION** Do not push forward or pull backward on the the active connector with the stent partially deployed. The active connector must be securely immobilized. Inadvertent movement of the the active connector may cause misalignment of the stent and possible damage.

#### **④ After Stent Deployment**

- Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
- Carefully remove the Electrocautery Stent Delivery System [HDS-10], guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion. (Place the inner sheath back into the outer sheath as the original state prior to removal.)
- Balloon dilatation inside the stent can be performed if judged necessary.

#### **12. Perform Routine Post-implant Procedures**

- Assess the good position of the stent and effective drainage. A stent may require up to 1 to 3 days to expand fully.
- Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- After implantation, the patient should remain on a soft diet until otherwise determined by the treating doctor.
- Observe the patient for development of any complications.

#### **13. Instructions for Removal of Niti-S Hot Nagi™ Stent (See Warnings)**

Grasp the retrieval string with forceps and/or collapse the proximal end of the Stent with snare then carefully and gently retrieve the Stent.

**Caution:** Do not apply excessive force to remove the stent as it may disconnect the retrieval string.

To reposition a Niti-S Hot Nagi™ Stent immediately after deployment, use forceps or a snare to grasp the retrieval string and gently adjust to the correct placement.

#### **14. Handling and Storage**

Store in a dry and cool place. Keep away from sunlight.

#### **Reuse Precaution Statement**

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your representative of Taewoong Medical Co., Ltd. For a single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Disposal Requirements:** The Electrocautery Stent Delivery System [HDS-10] of Niti-S Hot Nagi™ Stent must be properly sealed and disposed in compliance with the regulation of local or hospital at the end of its use.

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