Niti-S Esophageal Stent

Table of Contents

English	3
Français.....	9
Deutsch	15
Italiano	21
Español	27
Nederlands	33
Português	39
Türkçe.....	45
Čeština	51
Slovensky.....	56
Srpski	62
Dansk	67
Ελληνικά.....	73
Lietuvių.....	79
Polski	84
Română	90
Svenska	96
Norsk	101
Български.....	106

MRI Information

MR Conditional

Niti-S stent was determined to be MR-conditional.
Non-clinical testing demonstrated that the Niti-S Stent is MR Conditional according to ASTM F2503. A patient with this device can be scanned safely, immediately after placement under the following Conditions:

Static Magnetic Field

- Static magnetic field of 1.5-Tesla or 3-Tesla, only
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, Niti-S stent produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla / 64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris / 4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla / 128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	<u>1.5-Tesla</u>	<u>3-Tesla</u>
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.7-W/kg
Highest temperature change	+1.9°C	+2.4°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position Niti-S stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	7,695-mm ²	676-mm ²	7,984-mm ²	977-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

English

USER'S MANUAL

1. DESCRIPTION

The Niti-S Esophageal Stent consists of the implantable metallic stent and introducer system.

The stent is made of Nitinol wire. It is a flexible, fine mesh tubular prosthesis which has radiopaque markers on each end and at the center.

The device contains Nitinol, an alloy of nickel and titanium. Nickel may be released from the device. Nickel is classified as carcinogenic of category 2 and reproductive toxic of category 1B. Prior to implantation, patients should be counseled on the materials contained in the devices, as well as potential for carcinogenic and reproductive toxicity of these materials.


The Stent is loaded in introducer system and upon deployment the stent imparts an outward radial force on the luminal surface of the esophagus to establish patency.

The Anti-Reflex type Stent is designed to reduce or prevent reflux post-implantation.

Niti-S Full Covered Esophageal Stents used in benign stricture can be removed; (see Warnings).

Full Covered Esophageal Stents can be repositioned after deployment; (see Warnings).


A. Distal & Proximal Release Introducer

**Figure 1. Introducer System
(Distal release & Proximal release)**

- The introducer system accepts a 0.038" guidewire. The stent introducer system is passed over the guidewire into the esophagus.
- The stent is positioned appropriately using the X-ray markers for guidance under fluoroscopy.

B. TTS Introducer

**Figure 2. Introducer System
(TTS Introducer)**

- The TTS introducer system has a usable length of 160,180,220 cm
- TTS means **Through The endoScope**

2. PRINCIPLE OF OPERATION (DISTAL RELEASE & PROXIMAL RELEASE)

Distal release & TTS Introducer systems: Once the introducer is placed at its intended location, the outer sheath is pulled back by immobilizing the hub in one hand, grasping the Connector with the other hand, and gently sliding the Connector along the pusher towards the hub. Retraction of the outer sheath releases the stent. As soon as deployed, the stent expands to regain its original intended form and dimensions. While expanding, its expansion force pushes back the stricture in order to maintain esophageal luminal patency.

Proximal release systems: Once the introducer is placed at its intended location, the distal part of outer sheath is pushed forward by immobilizing the hub in one hand, grasping the Connector with the other hand, and gently sliding the Connector along the pusher towards the hub. Retraction of the outer sheath releases the stent. As soon as deployed, the stent expands to regain its original intended form and dimensions. While expanding, its expansion force pushes back the stricture in order to maintain esophageal luminal patency.

3. PATIENT POPULATION

Population generally over 18 years old, but decision belongs to the physicians

4. INTENDED USER

The device should only be used by trained and experienced physicians.

5. INTENDED PURPOSE

Maintaining esophageal luminal patency

6. INDICATION FOR USE

Indication for Use
Niti-S Esophageal Covered Stent [Both Bare-Type] is indicated for use in malignant strictures.

Figure 3. Stent Model

7. CLINICAL BENEFIT

The expected clinical benefit of the Niti-S Esophageal Stent is to maintain the lumen patency of the esophagus.

8. WARRANTY

Taewoong Medical Co., LTD. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Taewoong’s control directly affect the instrument and the results obtained from its use. Taewoong’s obligation under this warranty is limited to the repair or replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

* The user and/or patient should report to the competent authority of the Member State in which the user and/or patient is established and the manufacturer, if any serious incident that has occurred in relation to the device.

9. CONTRAINDICATION

The Niti-S Esophageal Stent is contraindicated for, but is not limited to:

- Patient with bleeding disorder.
- Strictures that do not allow passage of a guidewire.
- Any use other than those specifically outlined under indications for use.
- Removal or repositioning of fully deployed uncovered/bare Stents is contraindicated. (see Warnings).
- Suspected or impending perforation.

10. WARNINGS

- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the introducer system to organic solvent (e.g. alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- Niti-S Full Covered Stent cannot be removed when there is tumor in-growth/overgrowth/occlusion of the Stent lumen.
- Full Covered Stent may be repositioned immediately after deployment; see 18. Instructions for Removal of Full Covered Stents.
- Uncovered/bare Stents should not be removed once fully deployed; see Contraindications.
- Do not attempt to recapture/reload a stent once its deployment is advanced.
- Fully covered stents may be removed within 8 weeks. Stent removal shall be performed by a doctor according to the etiology of the benign stricture and the patient’s conditions.
- The risk of perforation and erosion into adjacent vascular structures or aorto-esophageal and arterio-esophageal fistulas may be increased with pre- or post-operative chemotherapy and radiation, longer implantation times, aberrant anatomy, and/or mediastinal contamination or inflammation.
- Silicone fully covered stents cannot be recaptured if the connector has been pulled beyond the pusher’s marker. Recapturing the stent in tortuous anatomy may damage the device. Recapturing more than once may also cause damages to the silicone membrane and/or the stent wire.

11. POTENTIAL COMPLICATIONS

Potential complications associated with the use and/or removal of Niti-S Stent may include, but are not limited to:

- Bleeding
- Inadequate expansion
- Pain
- Death (Other than due to normal disease progression)
- Deployment failure
- Perforation
- Stent misplacement
- Migration
- Stent occlusion
- Tumor overgrowth
- Tumor ingrowth
- Fever
- Foreign body sensation
- Sepsis
- Acute angulations
- Pneumonias
- Haematemesis
- Airway Compressions
- Reflux
- Food bolus impaction (lavage and debridement may be necessary on a periodic basis)
- Esophagitis
- Dysphagia
- Ulcerations
- Aspirations
- Stent fracture
- Mucosal tear
- Unsuccessful first removal attempt
- Esophageal avulsion
- Stridor requiring endotracheal intubation
- Fistula formation
- Esophagorespiratory fistula
- Impossibility to remove the stent
- Dislocation in stomach
- Cover breakdown with ingrowth of the mucosa
- Aorto and arterioesophageal fistula
- Erosion or perforation of stent into adjacent vascular structures
- Bilateral vocal cord paralysis

12. SSCP AND RESIDUAL RISK

Information on SSCP (Summary of safety and clinical performance), clinical benefit, Basic UDI and residual risk can be found on the website below.

Url: <https://ec.europa.eu/tools/eudamed>

Url: www.taewoongmedical.com

13. EQUIPMENT REQUIRED

- Distal or Proximal Release Introducer
 - Fluoroscope and/or endoscope
 - 0.038"/ 0.97 mm guidewire
 - Introducer sheath appropriately sized for stent and introducer system
- TTS Introducer
 - 0.035" (0.89mm) guidewire (preferably jag wire)
 - Introducer sheath appropriately sized for stent and introducer system
 - Endoscope system appropriately sized for instrument channel (8Fr or larger uncovered and covered, 3.7mm working channel)

14. PRECAUTIONS

Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- Care should be taken when removing the introducer system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging and the device should be inspected prior to use.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The Niti-S Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The Niti-S Stent is intended for single use only. Do not resterilize and/or reuse the device.

15. INSTRUCTIONS IN THE EVENT OF DAMAGE

WARNING: Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

16. PROCEDURE

① Examine stricture fluoroscopically and/or endoscopically.

- Carefully examine both the proximal and distal segment of stricture fluoroscopically.
- The Internal luminal diameter should be measured exactly with fluoroscope.

② Stent Size Determination

- Measure the length of the target stricture.
- Select a Stent size that is 20 to 40mm longer than the measured length of the stricture in order to cover fully both ends of lesion.
- Measure the diameter of the reference stricture - it is necessary to select a Stent which has an unconstrained diameter about 1 to 4mm larger than the largest reference target diameter, to achieve secure placement.

③ Stent Deployment Preparation

- The Niti-S Stent can be placed with the aid of fluoroscopy and/or endoscopy.
 - Pass a 0.038" (0.97mm) guidewire to the level of the stricture.
 - Pass a 0.035" (0.89 mm) guidewire to the level of the stricture.
- Under the fluoroscopy guidance, insert a guide wire across the stricture to where the stent introducer system will be placed over the guide wire.
 - Remove the stylet from the distal end of the introducer.
 - Ensure that the valve of connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.
 - Flush the inner lumen of introducer system.

④ Stent Deployment Procedure


Figure 4

PRECAUTION: Do not twist introducer system or employ a boring motion during the deployment as this may affect positioning and ultimate function of stent

A. Distal Release & TTS Introducer System

- Under the fluoroscopy and/or endoscopy guidance, position the introducer system to the center of target stricture exactly.
- Once the introducer system is in the correct position for deployment, unlock the proximal valve of the connector by turning the valve more than twice in an anti-clockwise direction. The stent is now ready for deployment
- To begin stent deployment, immobilize the hub in one hand and grasp the connector with the other hand. Gently slide the connector back along the pusher towards the hub.
- When the center X-ray marker reaches the center of target stricture, continue pulling back on the connector until the stent is fully deployed. (See figure 4, 5)


Figure 5

CAUTION Do not push forward or pull backward on the hub with the stent partially deployed. The hub must be securely immobilized. Inadvertent movement of the hub may cause misalignment of the stent and possible damage to the target or stricture.

B. Proximal Release System

- Under the fluoroscope and/or endoscopy guidance, position the introducer system to the center of target stricture exactly.
- Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Connector by turning the valve more than twice in an anti-clockwise direction. The stent is now ready for deployment
- To begin stent deployment, immobilize the Connector in one hand and grasp the hub with the other hand. Gently slide the hub forward along the 2nd inner catheter towards the Connector.
- When the center X-ray marker reaches the center of target stricture, continue forwarding toward the connector until the stent is fully deployed. (See figure 4, 6)


Figure 6

CAUTION Do not push forward or pull backward on the Connector with the stent partially deployed. The Connector must be securely immobilized. Inadvertent movement of the Connector may cause misalignment of the stent and possible damage to the esophagus.

⑤ After Stent Deployment

- Examine stent fluoroscopically to confirm expansion.
- Carefully remove the introducer system and the guidewire from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion. (Place the inner sheath back into the outer sheath as the original state prior to removal.)
- Balloon dilatation inside the Stent can be performed if the physician deems necessary.

17. PERFORM ROUTINE POST IMPLANT PROCEDURES.

- Assess the size and stricture of the Stent lumen. A Stent may require up to 1 to 3 days to expand fully.
- The Doctor should realize their experience and discretion in order to determine the appropriate drug regimen for each patient.
- After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- Observe the patient for development of any complications.

18. INSTRUCTIONS FOR REMOVAL OF NITI-S FULL COVERED STENTS (SEE WARNINGS)

Visually examine the Stent for any tumor in-growth/over-growth into the Stent lumen or whether the Stent is occluded. If the Stent lumen is clear, carefully remove using a forcep and/or snare. Grasp the retrieval string and/or collapse the proximal end of the Stent then carefully retrieve the Stent. If the Stent cannot be easily withdrawn, do not remove the Stent.

Caution: Do not allow excessive force to remove the stent as it may cause disconnect to the retrieval string.

To reposition a Niti-S Full Covered Stent immediately after deployment, use forceps or a snare to grasp the retrieval string and gently adjust to the correct placement.

Please note: the stent can only be repositioned and/or removed proximally.

19. REUSE PRECAUTION STATEMENT

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your Taewoong Medical Co., Ltd. representative. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and /or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

20. HANDLING AND STORAGE:

Store in a dry and cool place. Keep away from sunlight.

21. EXPECTED LIFETIME:

The device has been designed to function in the body for up to 1 year.

22. DISPOSAL REQUIREMENTS:

There are no special disposal requirements. After use, the introducer system and stent (when removed) must be properly sealed and disposed of in accordance with hospital, administrative and/or local government policy.

23. COMPOSITION

- Nitinol
- Platinum/Iridium
- Stainless Steel
- Silicone (applied only to covered stent)
- Polyester (applied only to covered stent)
- Nylon (applied only to Double & (Double)Anti-Reflux & Beta-2 stent)
- Polytetrafluoroethylene (applied only to (Double)Anti-Reflux & Beta-2 stent)

The Nitinol material is an alloy composed of nickel and titanium. Nickel may be released from the device and is classified as carcinogenicity of category 2 and reproductive toxicity of category 1B. Nickel patient exposure has been measured as 3.1640 µg/stent with a margin of safety of 6.95. Please note that a margin of safety greater than 1 for the potential leaching materials means that the medical device used under normal clinical conditions does not pose any toxicological risk to the adult patient population.

Besides nickel, the stent does not contain any other Carcinogenic, Mutagenic, Reprotoxic (CMR) and/or endocrine-disrupting substances in a concentration above 0,1 % weight by weight (w/w).

Prior to implantation, patients should be counseled on the materials contained in the devices, as well as potential for carcinogenic and reproductive toxicity of these materials.

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