


# Niti-S Biliary Covered Stent [Kaffes]

for benign and/or malignant biliary strictures



## Features

- **Long platinum radiopaque retrieval string:** The long platinum string helps easy removal
- **Radiopaque marker:** 3 (three) at both ends & 2 (two) in the middle

# Ordering Information

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TS0604F2	6	4	8.5	50
TS0605F2		5		
TS0606F2		6		
TS0607F2		7		
TS0608F2		8		
TS0804F2	8	4	8.5	50
TS0805F2		5		
TS0806F2		6		
TS0807F2		7		
TS0808F2		8		
TS1004F2	10	4	9	50
TS1005F2		5		
TS1006F2		6		
TS1007F2		7		
TS1008F2		8		

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BS0604F2	6	4	8.5	180
BS0605F2		5		
BS0606F2		6		
BS0607F2		7		
BS0608F2		8		
BS0804F2	8	4	8.5	180
BS0805F2		5		
BS0806F2		6		
BS0807F2		7		
BS0808F2		8		
BS1004F2	10	4	9	180
BS1005F2		5		
BS1006F2		6		
BS1007F2		7		
BS1008F2		8		